



Easing complex breathing

Philips Respironics BiPAP autoSV Advanced System One

BiPAP autoSV Advanced System One is our next step in the management of complex breathing patterns such as Cheyne-Stoke Respiration, CompSAS, and other types of Central Sleep Apnea.

BiPAP autoSV Advanced System One blends all proven features from the previous BiPAP AutoSV Advanced algorithm within the System One platform.

With auto EPAP, advanced back-up rate and sophisticated servo ventilation, the BiPAP autoSV Advanced System One may now provide a quieter and more comfortable therapy than before.

In addition to the EncorePro 2 complete reports, BiPAP autoSV Advanced System One provides physicians and caregivers with the possibility of an unequalled level of connectivity.

Our oximetry module and wireless modem work together to help clinicians create better patient management workflow.

PHILIPS

RESPIRONICS

sense and simplicity

Product specifications	
Dimensions	7" L x 5.5" W x 4" H 18 cm L x 14 cm W x 10 cm H
Weight	3 lbs; 1.36 kg (without humidifier)
Pressure range IPAP	4 cm H ₂ O to 25 cm H ₂ O
Pressure range EPAP	4 cm H ₂ O to 20 cm H ₂ O
Breath rate	Auto, 4-30bpm, off
Filters	Foam and optional ultra-fine
Device set-up	LCD/control wheel/push button
Data storage capacity (minimum)	Display: 7- and 30-day averages SD Card: >1 year On-board: up to 6 months
Compliance meter	Breathing detection
Altitude compensation	Automatic
Electrical requirements	100 – 240 VAC, 50/60 Hz
Humidification	Integrates with System One heated humidifier and pass-over humidifier
Warranty	Two years; replace with new

Ordering information		
Product description	Region	Part number
BiPAP autoSV Advanced with SD Card (device only)	International	IN951S
	France	FR955S
	Germany	DE952S
	LatAm	LA951S
	China	CN957S
BiPAP autoSV Advanced with SD Card (core package; includes device and heated humidifier)	Japan	JP956S
	International	IN951HS
	France	FR955HS
	Germany	DE952HS
	LatAm	LA951HS
	China	CN957HS
	Japan	JP956HS
Accessories		Part number
System One performance tube 15 (15 mm)		1074400
System One performance tube 22 (22 mm)		1032907
Device side panel		1063784
Water chamber		1066737
Power supply		1058190L
Carrying case		1063857
Foam filter kit		1063091
Ultra-fine filter kit (2 per pack)		1035442
System One heated humidifier		1056210
Link module (use with Sleep Therapy cable and System One devices)		1061644A
Power cord (5ft, 1.5m)	Region	Part number
Power cord	Australia	1039012
	Europe	1039014
	India	1060017
	Japan	1039027
	UK	1039026
Encore and SD Card		Part number
SD card (10 pack)		1063859
SD card mailer (10 pack)		1065146
SD card cover (attached to device)		1063858
SD card reader		1073915

Respirionics, autoSV, BiPAP, and Encore are trademarks of Respirionics, Inc. and its affiliates.

Please visit www.philips.com/respirionics



© 2011 Koninklijke Philips Electronics N.V. All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

CAUTION: US federal law restricts these devices to sale by or on the order of a physician.

Broudy NA 1/31/11 MCI 4104033 PN 1083045

Philips Healthcare is part of Royal Philips Electronics

Philips Respirionics International Headquarters
+33 1 47 28 30 82
Philips Respirionics Asia Pacific
+65 6882 5282
Philips Respirionics Australia
+61 (2) 9666 4444
Philips Respirionics United Kingdom
+44 800 1300 845

Philips Respirionics
1010 Murry Ridge Lane
Murrysville, PA 15668

Customer Service
+1 724 387 4000
800 345 6443 (toll free, US only)
www.philips.com/respirionics